

PROPRIETARY AND CONFIDENTIAL

THE INFORMATION CONTAINED IN THIS DRAWING IS THE SOLE PROPERTY OF TUTTNAUER LTD. ANY REPRODUCTION IN PART OR AS A WHOLE WITHOUT THE WRITTEN PERMISSION OF TUTTNAUER LTD. IS PROHIBITED.

125 ml



PREP. BY

DESCRIPTION




DATE

REV.

ECO






DETAIL B
SCALE 1 : 4



SECTION A-A

ITEM NO.	PART NUMBER	CATALOG NUMBER	125 ml/Q TY.
1	AJ0T11-0010	CHM507-0002	1
2	AJ4T52-0010-01	CMT507-0089	1
3	AM0T54-0030-02	BSK507-0003	3
4	DURAN- YOUTILITY bottle,with cap screw	Schott Duran 125 mml	147

	NAME	SIGNATURE	DATE	MODEL : 5075 ELV-MD		 Tuttnauer <small>Your Sterilization & Infection Control Partners</small>	A4 / ISO	 	
DESIGN	Yosi M		18.02.2014	REV NO.: 00	SCALE: 1:10				
DRAWN	Yosi M		18.02.2014	QUANTITY/UNIT: 1	QUANTITY/JOB:				DESCRIPTION: 5075 ELV with Youtility Scott Duran
APPROVED									
Q.A									
SHEET 1 OF 4			JOB NO.:	CAT. No.:		DRAWING NO.:	MATERIAL:		
				5075 ELV with Youtility Scott Duran					
Drawing folder : M:\Yosim\SolidWorks Work Space\5075 ELV with Youtility Scott Duran							WEIGHT:	218.26 Kg	

NOTES:
1) UNLESS OTHERWISE SPECIFIED: DIMENSIONS ARE IN MILLIMETERS.
2) DO NOT SCALE DRAWING.

PROPRIETARY AND CONFIDENTIAL

THE INFORMATION CONTAINED IN THIS DRAWING IS THE SOLE PROPERTY OF TUTTNAUER LTD. ANY REPRODUCTION IN PART OR AS A WHOLE WITHOUT THE WRITTEN PERMISSION OF TUTTNAUER LTD. IS PROHIBITED.

Prep.by


Description

Date





Rev. No.

ECO

250 ml






DETAIL D
SCALE 1 : 6



SECTION C-C

ITEM NO.	PART NUMBER	CATALOG NUMBER	250 ml/QTY.
1	AJ0T11-0010	CHM507-0002	1
2	AJ4T52-0010-01	CMT507-0089	1
3	AM0T54-0030-02	BSK507-0003	3
4	DURAN- YOUTILITY bottle,with cap screw	Schott Duran 250 mm	111

	NAME	SIGNATURE	DATE	MODEL : 5075 ELV-MD		 Tuttnauer Your Sterilization & Infection Control Partners	A4 / ISO	 
DESIGN	Yosi M		18.02.2014	REV NO.: 00	SCALE: 1:12			
DRAWN	Yosi M		18.02.2014	QUANTITY/UNIT: 1	QUANTITY/JOB:	DESCRIPTION: 5075 ELV with Youtility Scott Duran		
APPROVED								
Q.A								
SHEET 2 OF 4			JOB NO.:	CAT. No.:		DRAWING NO.:	MATERIAL:	
						5075 ELV with Youtility Scott Duran		
Drawing folder : M:\Yosim\SolidWorks Work Space\5075 ELV with Youtility Scott Duran							WEIGHT:	218.26 Kg

NOTES:

1) UNLESS OTHERWISE SPECIFIED: DIMENSIONS ARE IN MILLIMETERS.
2) DO NOT SCALE DRAWING.

PROPRIETARY AND CONFIDENTIAL

THE INFORMATION CONTAINED IN THIS DRAWING IS THE SOLE PROPERTY OF TUTTNAUER LTD. ANY REPRODUCTION IN PART OR AS A WHOLE WITHOUT THE WRITTEN PERMISSION OF TUTTNAUER LTD. IS PROHIBITED.

Prep.by





Description

Date

Rev. No.

ECO


500 ml



DETAIL F
SCALE 1 : 6

SECTION E-E

ITEM NO.	PART NUMBER	CATALOG NUMBER	500 ml/QTY.
1	AJ0T11-0010	CHM507-0002	1
2	AJ4T52-0010-01	CMT507-0089	1
3	AM0T54-0030-02	BSK507-0003	3
4	DURAN- YOUTILITY bottle,with cap screw	Schott Duran 500 mml	69

	NAME	SIGNATURE	DATE	MODEL : 5075 ELV-MD	 A4 / ISO	
DESIGN	Yosi M		18.02.2014	REV NO.: 00		
DRAWN	Yosi M		18.02.2014	SCALE: 1:12	5075 ELV with Youtility Scott Duran	
APPROVED				QUANTITY/UNIT: 1		
Q.A				QUANTITY/JOB:	5075 ELV with Youtility Scott Duran	
SHEET 3 OF 4			JOB NO.:	CAT. No.:	5075 ELV with Youtility Scott Duran	
Drawing folder : M:\Yosim\SolidWorks Work Space\5075 ELV with Youtility Scott Duran					WEIGHT:	218.26 Kg

NOTES:

1) UNLESS OTHERWISE SPECIFIED: DIMENSIONS ARE IN MILLIMETERS.
2) DO NOT SCALE DRAWING.

PROPRIETARY AND CONFIDENTIAL

THE INFORMATION CONTAINED IN THIS DRAWING IS THE SOLE PROPERTY OF TUTTNAUER LTD. ANY REPRODUCTION IN PART OR AS A WHOLE WITHOUT THE WRITTEN PERMISSION OF TUTTNAUER LTD. IS PROHIBITED.

Prep.by

Description

Date

Rev. No.

ECO

-





See Sheet1

-


-

-




1000 ml



SECTION G-G



ITEM NO.	PART NUMBER	CATALOG NUMBER	1000 ml/QTY.
1	AJ0T11-0010	CHM507-0002	1
2	AJ4T52-0010-01	CMT507-0089	1
3	DURAN- YOUTILITY bottle,with cap screw	Schott Duran 1000 mml	38
4	AM0T54-0010-02	BSK507-0001	2

	NAME	SIGNATURE	DATE	MODEL : 5075 ELV-MD		 Tuttnauer <small>Your Sterilization & Infection Control Partners</small>	A4 / ISO	 
DESIGN	Yosi M		18.02.2014	REV NO.: 00	SCALE: 1:12			
DRAWN	Yosi M		18.02.2014					
APPROVED								
Q.A				QUANTITY/UNIT: 1	QUANTITY/JOB:	DESCRIPTION: 5075 ELV with Youtility Scott Duran		
SHEET 4 OF 4			JOB NO.:	CAT. No.:		DRAWING NO.:		MATERIAL:
5075 ELV with Youtility Scott Duran								
Drawing folder : M:\Yosim\SolidWorks Work Space\5075 ELV with Youtility Scott Duran							WEIGHT:	218.26 Kg

NOTES:

1) UNLESS OTHERWISE SPECIFIED: DIMENSIONS ARE IN MILLIMETERS.
2) DO NOT SCALE DRAWING.