

**PROPRIETARY AND CONFIDENTIAL**

THE INFORMATION CONTAINED IN THIS DRAWING IS THE SOLE PROPERTY OF TUTTNAUER LTD. ANY REPRODUCTION IN PART OR AS A WHOLE WITHOUT THE WRITTEN PERMISSION OF TUTTNAUER LTD. IS PROHIBITED.

PREP. BY



DESCRIPTION

DATE




REV.

ECO



125 ml



SECTION A-A



ITEM NO.	PART NUMBER	DESCRIPTION	CATALOG NUMBER	125 ml/Q TY.
1	AG0T51-0010	Tray Holder	TRH507-0002	1
2	DURAN- YOUTILITY bottle,with cap screw	Schott Duran 125 mml	Schott Duran 125 mml	154

	NAME	SIGNATURE	DATE	MODEL :	5075 shelf	 A4 / ISO	
DESIGN	Yosi M		17.02.2014	REV NO.:	00		
DRAWN	Yosi M		17.02.2014	SCALE:	1:10		
APPROVED				QUANTITY/UNIT:	1		
Q.A				QUANTITY/JOB:		DESCRIPTION:	
SHEET 1 OF 4			JOB NO.:	CAT. No.:		DRAWING NO.:	MATERIAL:
5075 Loading shelf with Youtility Scott Duran							
Drawing folder : M:\Yosim\SolidWorks Work Space\5075 Loading shelf with Youtility Scott Duran							WEIGHT: Kg

**NOTES:**

1) UNLESS OTHERWISE SPECIFIED: DIMENSIONS ARE IN MILLIMETERS.  
2) DO NOT SCALE DRAWING.



PROPRIETARY AND CONFIDENTIAL

THE INFORMATION CONTAINED IN THIS DRAWING IS THE SOLE PROPERTY OF TUTTNAUER LTD. ANY REPRODUCTION IN PART OR AS A WHOLE WITHOUT THE WRITTEN PERMISSION OF TUTTNAUER LTD. IS PROHIBITED.

Prep.by

Description

Date

Rev. No.

ECO

-




See Sheet1

-



-

-




250 ml



SECTION C-C



ITEM NO.	PART NUMBER	DESCRIPTION	CATALOG NUMBER	250 ml/QTY.
1	AG0T51-0010	Tray Holder	TRH507-0002	1
2	DURAN- YOUTILITY bottle, with cap screw	Schott Duran 250 mml	Schott Duran 250 mml	91

	NAME	SIGNATURE	DATE	MODEL :	5075 shelf		 <b>Tuttnauer</b> <small>Your Sterilization &amp; Infection Control Partners</small>	A4 / ISO	 	
DESIGN	Yosi M		17.02.2014	REV NO.:	00	SCALE:				1:10
DRAWN	Yosi M		17.02.2014	QUANTITY/UNIT:	QUANTITY/JOB:					DESCRIPTION:
APPROVED				1						
Q.A										
SHEET 2 OF 4			JOB NO.:	CAT. No.:		DRAWING NO.:		MATERIAL:		
				5075 Loading shelf with Youtility Scott Duran						
Drawing folder : M:\Yosim\SolidWorks Work Space\5075 Loading shelf with Youtility Scott Duran								WEIGHT:	Kg	

NOTES:



1) UNLESS OTHERWISE SPECIFIED: DIMENSIONS ARE IN MILLIMETERS.  
2) DO NOT SCALE DRAWING.





**PROPRIETARY AND CONFIDENTIAL**

THE INFORMATION CONTAINED IN THIS DRAWING IS THE SOLE PROPERTY OF TUTTNAUER LTD. ANY REPRODUCTION IN PART OR AS A WHOLE WITHOUT THE WRITTEN PERMISSION OF TUTTNAUER LTD. IS PROHIBITED.


**500 ml**




SECTION A-A



need to remove



ITEM NO.	PART NUMBER	DESCRIPTION	CATALOG NUMBER	500 ml/QTY.
1	AG0T51-0010	Tray Holder	TRH507-0002	1
2	DURAN- YOUTILITY bottle,with cap screw	Schott Duran 500 mml	Schott Duran 500 mml	30

	NAME	SIGNATURE	DATE	MODEL :	5075 shelf	 A4 / ISO	
DESIGN	Yosi M		17.02.2014	REV NO.:	00		
DRAWN	Yosi M		17.02.2014	SCALE:	1:10	DESCRIPTION:	
APPROVED				QUANTITY/UNIT:	1		
Q.A				QUANTITY/JOB:		DRAWING NO.:	
SHEET 3 OF 4			JOB NO.:	CAT. No.:			
5075 Loading shelf with Youtility Scott Duran						MATERIAL:	
Drawing folder : M:\Yosim\SolidWorks Work Space\5075 Loading shelf with Youtility Scott Duran						WEIGHT: Kg	

**NOTES:**

1) UNLESS OTHERWISE SPECIFIED: DIMENSIONS ARE IN MILLIMETERS.  
2) DO NOT SCALE DRAWING.



**PROPRIETARY AND CONFIDENTIAL**

THE INFORMATION CONTAINED IN THIS DRAWING IS THE SOLE PROPERTY OF TUTTNAUER LTD. ANY REPRODUCTION IN PART OR AS A WHOLE WITHOUT THE WRITTEN PERMISSION OF TUTTNAUER LTD. IS PROHIBITED.

Prep.by

Description

Date

Rev. No.

ECO


-

See Sheet1

-




-

-





1000 ml

need to remove



SECTION D-D

ITEM NO.	PART NUMBER	DESCRIPTION	CATALOG NUMBER	1000 ml/QTY.
1	AG0T51-0010	Tray Holder	TRH507-0002	1
2	DURAN- YOUTILITY bottle,with cap screw	Schott Duran 1000 mml	Schott Duran 1000 mml	21

	NAME	SIGNATURE	DATE	MODEL :	5075 shelf	 A4 / ISO	
DESIGN	Yosi M		17.02.2014	REV NO.:	00		
DRAWN	Yosi M		17.02.2014	SCALE:	1:10		
APPROVED				QUANTITY/UNIT:	1		
Q.A				QUANTITY/JOB:		DESCRIPTION:	
SHEET 4 OF 4			JOB NO.:	CAT. No.:	5075 Loading shelf with Youtility Scott Duran		
Drawing folder : M:\Yosim\SolidWorks Work Space\5075 Loading shelf with Youtility Scott Duran							WEIGHT: Kg

**NOTES:**

1) UNLESS OTHERWISE SPECIFIED: DIMENSIONS ARE IN MILLIMETERS.  
2) DO NOT SCALE DRAWING.