




Model: 2840ELC-G-D
Including steam generator



Front view



Left side view




Back view

Model	Chamber diameter x depth	Volume
2840ELC-G-D	280x 400	28.5 Liter
2840ELC-PVG-D	280x 400	28.5 Liter

Electrical option	Standard
Voltage	208V/230V
Electric Power	2.2kW
Amperage	10A
Phase	1 Ph

No.	Description	Size	Minimum Pressure
A	Over flow	1/2" -Male	-
B	Tap water inlet	1/2" -Male	-
C	Drain	1/2" -Male	-
D	Air pressure inlet	1/2" -Male	3 Bar
E	Power supply	-	-
F	Ethernet port	RJ-45	-
G	Mineral water inlet	-	-



Top view

Installation drawing