




Model: 2840 EL-D



Front View




Left side View



Back View


Model: 2840 EL-PV-D (Includ vavuum pump)



Back View

Model	Chamber diameter x depth	Volume
2840 EL-D	280x 400	28.5 Liter

Electrical specification	Standard
Voltage	208V/230V
Electric Power	2.2kW
Amperage	10A
Phase	1 Ph



Top View

Installation drawing

No.	Description	Size	Minimum Pressure
A	Over Flow (PV model only)	1/2" -Male	-
B	Mineral water inlet	-	-
C	Power supply	-	-
D	Ethernet port	RJ-45	-